BENTOGRAN

Granular bentonite

ALLERGEN	GMO
FREE.	FREE

🔽 🗊 💓 🐙

→ TECHNICAL DESCRIPTION

Bentogran is a bentonite of easy utilization because it dissolves promptly. The granular formulation enables to use it without creating any dust, with great advantage for the personnel in charge of dissolving it.

Bentogran suspensions are ready to be used after 2 hours, unlike those of traditional bentonites which require 8-12 hours of swelling.

It absorbs a high quantity of water, with a swelling index that highlights its exceptional adsorbent and super deproteinizing qualities, particularly appreciated in the production of quality wines.

-> COMPOSITION AND TECHNICAL CHARACTERISTICS

Activated bentonite.

→ DOSAGE

From 10 to 60 g/hL.

→ INSTRUCTIONS FOR USE

Rehydration in 10 parts of cold water will produce a creamy gel.

→ STORAGE AND PACKAGING

Store in a cool dry place, away from direct sunlight and heat.

1 kg net packs in cartons containing 20 kg. 25 kg net bags.

